

BELT GRINDER

Several Purpose Belt Grinder



For flat, coarse grinding prior to fine grinding. The specimens to be ground in 60,80 or 120 grit after cutting operation or after mounting in mould to remove surface imperfections. Useful step to remove light metallurgical damage in the specimen surface induced during cut-off operation.



- Sturdy Table Top Model for Pre-Grinding of Metallurgical Specimen
- ½ HP Single Phase Motor , 1440 RPM
- 100mm x 915mm Endless Belt
- Constant Tension Spring Loaded mechanical system
- Aluminium Idler Pulley Mechanism
- Simplified Belt release mechanism
- Easy removable Dust cover system
- Size: 660 x 430 x 290 mm

Optional Spares

On-Off Switch

Bainline BG, Heavy-Duty belt grinder,
Twin Belt Grinder and Heavy-Duty Spectro grinder
also available



Chennai Metco Pvt Ltd

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India.
Ph: 044-2688 0737, 044-2688 0005 Fax: (044) 2688 0850,
e-mail: mktg@chennaietco.com, web : www.chennaietco.com