



BAINCUT - L

Large Abrasive Cutter



A Sturdy floor model general purpose cutter. Large space in the cutting zone offers flexibility to use wide variety of fixtures. The longitudinal table movement across cutting area (Y-movement), side movement (X-movement) offers extra convenience. Suitable for variety of applications including gears, lengthy components.

- Sturdy floor model cutter for Metallographic sectioning
- Rust proof Steel construction
- Cutting Capacity upto 120 mm dia or 100 x 150 mm
- 7.5 HP, 3 Phase Motor
- Spindle speed 2880 rpm
- cut-off wheel dia 12"/14"
- Splash proof, corrosion resistant see-through hood
- Cooling by Four high flow water jets to provide optimum cooling
- Coolant Tank with 190 ltrs capacity
- 65 mm drain for coolant recirculation
- 1/3 HP Coolant Recirculation Pump
- Flash cleaning attachment with hook for mounting
- Fluorescent light in the working area provide clear illumination
- Start button with key for safety
- T Slot Bed
- Cam vise
- Y Movement (Longitudinal sectioning)
- Machine Size: 1020x1260x1030 mm
- Stand size : 1320x1270x810 mm



Automation also available

Chennai Metco Pvt Ltd

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India.
Ph: 044-2688 0737, 044-2688 0005, Fax: (044) 2688 0850,
e-mail: mktg@chennaietco.com, web : www.chennaietco.com