

Baincut L Auto

LARGE HEAVY-DUTY CUTTING M/C



Chennai Metco

www.chennaimetco.com

- Sturdy floor model cutter for Metallographic sectioning
- Rust proof Steel construction
- Cutting Capacity up to 120 mm dia or 100 x 150 mm
- 7.5 HP, 3 Phase Motor
- Spindle speed 2880 rpm
- Cut-off wheel dia 12"/14"
- Cutting action Z and Y movements
- Y Movement Automatic Cutting
- Electronic Automation with servo motor
- Splash proof, corrosion resistant doors and see-through hood
- Cooling by Four high flow water jets to provide optimum cooling
- Coolant Tank with 150 ltrs capacity
- One large 65 mm drain for coolant recirculation
- 1/2 HP Recirculation Pump
- Flash cleaning attachment with hook for mounting
- Fluorescent light in the working area provides clear illumination
- T-Slot Bed, Camvise
- Y Movement (Longitudinal sectioning)



Size

Machine Size: 1320 x 1350 x 810 mm

Stand size : 1320 x 1270 x 810 mm

Baincut Large Manual also available

Chennai Metco Pvt Ltd

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India.

Ph: 044-2688 0737, 04-2688 0005, Fax: (044) 2688 0850

e-mail: mktg@chennaietco.com, web : www.chennaietco.com