

Bainmount - Pneumatic

Pneumatic Hot Mounting Press

This machine is used in metallurgical lab for mounting the sample and for the purpose of viewing the micro structure, surface finish and edge retention.



BAINMOUNT PNEUMATIC SPECIFICATIONS

- Corrosion resistant steel Cabinet
- Pneumatic operation System
- Digital Temperature Indicator and digital timer
- Quick Ram return to speed up reloading
- 1500 Watts mould Heater
- Dovetail slot for easy change of Mould Assemblies
- Water Cooling System (manual)
- Timer-Buzzer indication for heating and cooling cycles
- Wide range of sample sizes can be accommodated
- Two moulds can be simultaneously prepared using Dual Mould Spacer
- Provision to accommodate 1 1/2" & 2" Mould Sizes
- Size : 580 x 630x 600 mm
- Weight - 78kg



Optional Accessories

- Mould Assembly 1 1/2"
- Mould Assembly 2"
- Water Cooling (Automatic)
- Recirculating Water Cooling system

Hydraulic also available

Chennai Metco Pvt Ltd

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India.
Ph: 044-2688 0737, 044-2688 0005 Fax: (044) 2688 0850,
e-mail: mktg@chennaietco.com, web : www.chennaietco.com