



BAINCUT- M Auto

Metallography

Abrasive Cutting Machine



This machine is used in the metallurgical lab for cutting the selected portion of the component as cut sample for mounting it in the mould and viewing the case depth microstructure without affecting its original property.



- Sturdy table model cutter for Metallography sectioning of Metals, Ceramic and mineral samples
- Rust proof Steel construction
- 3HP, 3 Phase Motor
- Cutting capacity upto 60 mm(Standard size)
- Spindle speed 2880 rpm with automatic cutting
- Cut-off wheel dia 10"
- Splash Proof, corrosion resistant with see through hood
- Cooling by two water jets to provide optimum cooling
- Inbuilt movable recirculation coolant tank with 50 ltrs capacity
- One large drain for coolant recirculation
- 1/3 HP Coolant Pump
- Start button with key for safety
- Small opening in the side wall for long jobs
- Fluorescent light in the working area provides Safe and clear illumination

Size

- Machine Size : 640 x 920 x 560 mm
- Stand Size : 850 x 760 x 730 mm

Baincut — Manual also available

Chennai Metco Pvt Ltd

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India.

Ph: 044-2688 0737,044-2688 0005 Fax: 044-2688 0850

e-mail: mktg@chennaietco.com, web : www.chennaietco.com