



BAINPOL ET

SINGLE DISC TABLE TOP GRINDER / POLISHER



Fully moulded FRP caver. Quick change disc design enables discs with different emery grits/polish cloths to be used in desired steps.



www.chennaietco.com

Specifications

- Rigid table top corrosion proof FRP cabinet
 - 8" / 10" / 12" Single Disc Polisher
- Single Speed (720 rpm) 0.5HP high torque AC motor
Corrosion resistant wash bowl, free flow drain
Flexible Water Jet with control valve
Aluminium Disc 12" Dia (Interchangeable)
Holding Band
S.S Holding Ring (Press Type)

Optional Accessories

- Magnetic substrate disc to attach to any conventional platen self adhesive backing
- A Set of Ferromagnetic carrier sheet for sand paper/cloths

Bainpol ETD , VTD (variable speed double disc),
VFD (variable speed floor model) ,
PM Auto (Polishing machine fully auto) and SA available

Chennai Metco Pvt Ltd

SP100A, Ambattur Industrial Estate,
Chennai 600 058.

Ph: 044-2688 0737/ 0005, Fax: 044-2688 0850
e-mail: mktg@chennaimetco.com

