

Bainpol ET


Single Disk Single Speed



For Grinding/Polishing of the Specimen

Single speed Polisher/Grinder with fully moulded FRP body. Quick change disc design enables discs with different emery grits/polish cloths to be used in desired steps.

- 8"/12" Single disc polisher
- Rigid table top Model with corrosion proof FRP Cabinet
- Single Speed 0.5HP high torque AC motor
- Corrosion resistant wash bowl, free flow drain
- Flexible Water Jet with control valve
- Size: 375 X 655 X 330 mm
- Aluminium Disc 8"/12" Dia (Interchangeable)
- Holding Spring
- S.S Holding Ring (Press Type)
- Size: 360 x 660 x 310mm



VTD,VT, ETD, VFD and SA also available



Chennai Metco Pvt Ltd

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India.
Ph: 044-2688 0737, 044-2688 0005, Fax: (044) 2688 0850,
e-mail: mktg@chennaietco.com, web : www.chennaietco.com