



# BAINPOL ETD

DOUBLE DISC TABLE TOP GRINDER / POLISHER



**Fully moulded FRP caver. Quick change disc design enables discs with different emery grits/polish cloths to be used in desired steps.**



[www.chennaietco.com](http://www.chennaietco.com)

## Specifications

- Rigid table top corrosion proof FRP cabinet
- Two Speed System ( 1410 to 1460 rpm and 700 to 720rpm)
- 8" / 10" / 12" double disc polisher
- 0.5HP AC Motor
- Corrosion resistant wash bowl
- Free flow drain system
- Flexible water jet with control valve
- Aluminium Disc 8" Dia (Interchangeable)
- S.S Holding Ring (Press Type)
- Holding Band

## Optional Accessories

- Magnetic substrate disc to attach to any conventional platen self adhesive backing
- A Set of Ferromagnetic carrier sheet for sand paper/cloths

Bainpol ET , VTD (variable speed double disc),  
VFD ( variable speed floor model) ,  
PM Auto (Polishing machine fully auto) and SA available

## Chennai Metco Pvt Ltd

SP100A, Ambattur Industrial Estate,  
Chennai 600 058.

Ph: 044-2688 0737/ 0005, Fax: 044-2688 0850  
e-mail: [mktg@chennaimetco.com](mailto:mktg@chennaimetco.com)

