

Bainmount - P Auto

Automatic Mounting Press Pneumatic

High speed heater with digital display. Suitable for labs with moderate number of samples. Automation offers reduced manual intervention and improved consistency.



Chennai Metco

- Steel Cabinet with built in Pneumatic System
- Digital Temperature indicator and Digital timer with LCD Display
- Quick Ram return to speed up reloading
- 1200 Watt mould Heater
- Water Cooling System
- Thermostatically controlled 1500 Watts Mould heater
- Buzzer indication for cooling cycles
- Wide Range of Sample size can be accommodated
- With Interchangeable 1¼" , 1½" 2" mould sizes
- Size : 580 x 630x 600 mm
- Weight - 78kg



Optional Accessories

- Mould Assembly 1½"
- Mould Assembly 2"

Hydraulic also available

Chennai Metco Pvt Ltd

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India.

Ph: 044-2688 0737, 044-2688 0005, Fax: (044) 2688 0850

e-mail: mktg@chennaietco.com, web : www.chennaietco.com