

METALLOGRAPHY

BAINPOL AUTO

Automatic polishing machine
practical and economic.

Can be the choice for the laboratory
looking for reproducible results.

Pneumatic assisted.

Can hold multiple moulds (Upto 5
moulds)



Chennai Metco

BAINPOL Auto

Automatic Polishing Machine

Practical and economic auto polishing machine. Can be the choice for the laboratory looking for reproducible results. Pneumatic assisted. Can hold multiple moulds (upto 5 moulds) - 1 1/4" dia (standard), 1 1/2" dia. Disc size 10" 912" also available)

- Variable Speed 50- 600 RPM
- 10" Single Aluminum Disc Hold up to 5 Moulds (1 1/4" standard)
- 1 HP High torque AC Motor
- 1/2 HP Independent Polishing Head Motor
- Polisher Head Variable Speed : 50 to 200 RPM
- Electronic Control with LC - Display with Digital timer
- Variable force : 1 to 4 bar
- With Interchangeable 1 1/4" , 1 1/2" 2" Mould sizes
- Single Pressure Specimen loading, Pneumatically applied
- Power Supply : 110V (Single Phase + Neutral + Earth)
- Size: 700 X 700 X 750 mm



Chennai Metco Pvt Ltd

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India. Ph: 044-2688 0737, 044-2688 0005, Fax: (044) 2688 0850, e-mail: mktg@chennaietco.com, www.chennaietco.com