


PORTABLE

MICROSCOPE

We make available portable Wide Stand Microscope series with coaxial incident light illumination. The illumination comes down from the inside of the microscope.

This illumination system is suitable for observing Microstructure onsite.



- Magnification : 100X
- Eye piece : WF10X
- Objective: 10X
- Illumination: Opac Illuminator w/small Bulb 4.8 V, 0.5A and 1 meter long cord
- Dry Cells : 4 cell (1.5V X 4 = 6C) in case continuous employment abt 5.5h

Chennai Metco

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India. Ph: 044-2688 0737, Fax: (044) 2688 0850, e-mail: mktg@chennaietco.com, web : www.chennaietco.com

www.chennaietco.com