


SEMI AUTOMATIC

# Grinding/Polishing Machine



**Chennai Metco Pvt Ltd.,**

[www.chennaimetco.com](http://www.chennaimetco.com)

# BAINPOL SA

## Semi Automatic Polishing Machine



Practical and economic auto polishing machine. Can be the choice for the laboratory looking for reproducible results. Can hold multiple moulds 1 ¼", 1 ½" and 2". Non complicated design keeping Indian context in mind.

### Description:

Holds upto 3 moulds 1¼" / 1½" / 2" Dia  
Rigid Table Top Corrosion Proof FRP Cabinet  
0.5 HP AC high torque Motor  
Imported Drive for smooth speed change  
Corrosion resistant wash bowl with free flow drain design  
Flexible Water Cooling System  
Interchangeable Disc Design  
Variable Speed 50 - 500 RPM  
12" Aluminium Alloy Disc(Anodized)  
Holding spring  
Size: 435 x 750 x 520 mm  
1 1/4" holder comes with Machine as a Standard

### Polishing Head

Independent Powered Polish Head  
Digital Timer  
Single Speed AC Motor  
(110 V, 40W, 100 RPM)



### Chennai Metco Pvt Ltd

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India.

Ph: 044-2688 0737, 044-2688 0005, Fax: (044) 2688 0850,

e-mail: [mktg@chennaietco.com](mailto:mktg@chennaietco.com), [www.chennaietco.com](http://www.chennaietco.com)