

# Bainline - VSG

## Variable speed spectro grinder

For perfect flat grinding  
for perfect spectroscopy  
results.

These specialised units  
generate appropriate  
surface for spectroscopy.

Vacuum suction and  
magnetic holder  
complete setup.

Variable speed option for  
multiple applications.



This machine is designed as per International Standards to produce the surface of quality suitable for subsequent spectroscopic analysis

**Chennai Metco**

- 14" Vacuum Disc & Emery Disk Holder
- Floor Model
- Disc Speed 100 to 2880 RPM
- 1 HP Power
- VACUUM SUCTION SYSTEM
- IMPORTED FREQUENCY DRIVE(For Variable Speed 100-2880rpm)
- Magnetic Holder
- Power Supply : 415 V / 50Hz ( R + Y + B + N + E )
- Size - 1050X650x1000mm
- 350x40mm Grinding Disc Grit Size 36



A Vacuum Suction dust collector is provided for accumulation of wastes

Bainline SSG  
(Single Speed)  
also available



## **Chennai Metco Pvt Ltd**

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India.  
Ph: 044-2688 0737, 044-2688 0005 Fax: 044 2688 0850,  
e-mail: [mktg@chennaietco.com](mailto:mktg@chennaietco.com), web : [www.chennaietco.com](http://www.chennaietco.com)