

# TWIN BELT GRINDER

# WET



## FOR FLAT COARSE GRINDING

For flat, coarse grinding prior to fine grinding. The specimens to be ground in 60,80 or 120 grit after cutting or moulding operation to remove surface imperfections. Useful step to remove light metallurgical damage in the specimen surface induced during cut-off operation.

[www.chennaietco.com](http://www.chennaietco.com)

- Grinding belts on both sides
- water cooling
- Mechanism for belt centering
- SS cover
- Motor: 1 HP
- Belt Speed: 1500 RPM
- Grinding belts: 100 x 920 mm, 100 X 915 mm
- Water Connection
- Drain : Right and Left, Ø 40 mm
- Water Valve: Right & Left
- Dimension.W x H x D 600 X 750 X 270



**Bainline GP, Bainline - H (Heavy-Duty)  
and Spectro Grinder Available**



## **Chennai Metco Pvt Ltd**

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India.

Ph: 044-2688 0737, 044-2688 0005 Fax: 044 - 2688 0850

e-mail: [mktg@chennaietco.com](mailto:mktg@chennaietco.com), web : [www.chennaietco.com](http://www.chennaietco.com)