



BAINPOL - VFD

Heavy duty floor model
Grinder/Polisher



Fully moulded FRP body. Quick change disc design enables discs with different emery grits/polish cloths to be used in desired steps. Particularly suitable for High volume/Large specimen preparation. Includes compartments for storing consumables.

- 12" double disc polisher
- Rigid floor model corrosion proof steel cabinet
- Variable Speed AC Drive (Imported)
- 0.5HP imported AC Motor (Individual Drive)
- Corrosion Resistant wash bowl with free
- Flow drain design
- Flexible water jet with control valve
- Aluminium Disc 12"Dia (Inter Changeable)
- Holding Spring
- Holding Ring (Press Type)
- Size: 870x750x1020mm



1. Digital Display
2. Disc Run
3. Disc Stop
4. Forward/Reverse
5. RPM Decrement
6. RPM Increment

Also Available (Table top model)

- VTD : Double Disc Variable Speed
- VT : Single Disc Variable Speed
- ETD : Double Disc Single Speed
- ET : Single Disc Single Speed
- SA : Semi Auto

Chennai Metco Pvt Ltd

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India.

Ph: 044-2688 0737, 044-2688 0005 Fax: (044) 2688 0850

e-mail: mktg@chennaietco.com, web : www.chennaietco.com