

# Bainpol VT

Single Disk Variable Speed



## For Grinding/Polishing of the Specimen

Variable speed Polisher/Grinder with fully moulded FRP body. Quick change disc design enables discs with different emery grits/polish cloths to be used in desired steps.

- 8"/12" Single disc polisher
- Variable speed from 50-1000 RPM
- Rigid table top Model, corrosion proof FRP cabinet
- 0.5HP High torque Motor
- Imported Drive for smooth speed change
- Corrosion Resistant wash bowl of free flow drain design
- Flexible water jet with control valve
- Aluminium Disc 8"/12" Dia (Interchangeable)
- Holding Spring
- Stainless Steel Holding Ring (Press Type)
- Size: 360 x 660 x 310mm



## Optional Feature

Front Control Panel with feather touch buttons

VTD, ET, ETD, VFD and SA also available



## **Chennai Metco Pvt Ltd**

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India.  
Ph: 044-2688 0737, 044-2688 0005, Fax: (044) 2688 0850,  
e-mail: [mktg@chennaietco.com](mailto:mktg@chennaietco.com), web : [www.chennaietco.com](http://www.chennaietco.com)