

BAINPOL - VT/VTD


Single/Double disc variable speed Polisher/grinder

Variable speed polisher/grinder fully moulded FRP body. Quick change disc design enables discs with different emery grits/polish cloths to be used in desired steps.

- 0.5HP high torque motor
- Imported AC drive
- 8" disc dia. standard
- 12" dia options
- LCD display, 50 to 1000 rpm
- Size VT :360 x 660 x 310 mm
- VTD : 710 x 660 x310 mm
- Detachable coolant Jet



PMV008



PMV009



Feather Touch Control Panel

- | | |
|-------------------------|--------------------|
| a. RPM/Time Increment | d. Disc Stop/Run |
| b. RPM/Time Decrement | e. Power On/Off |
| c. Disc Forward/Reverse | f. Digital Display |

Single speed polisher/grinder also available