


BAINCUT

XL


AUTO



Automatic bigger size fully loaded cutting machine with top closed design. Voluminous internal space to keep big components.

www.chennai metco.com

Sturdy floor model cutter for sectioning
 Rust proof steel construction
 Cutting capacity upto 150mm or 100 X 150 mm
 10 HP, 3 Phase Motor
 Spindle speed 1440 RPM
 Cut-off wheel dia 16"/20"
 Cutting action automatic - 'Z' Movement
 The machine also offers manual mode
 Splash proof, corrosion resistant doors
 Cooling by multiple high flow water jets to provide optimum cooling
 Coolant Tank with 250 ltrs capacity
 ½ HP Recirculating pumps
 Two large drains for draining the coolant
 Flash cleaning attachment with hook for mounting
 Fluorescent light in the working area provides safe & clear illumination
 Start button with key for safety
 T Slot bed, Cam Vise, X Movement, Y Movement
 Machine Size : 2100 X 1750 X 2070 mm



Optional Accessories

Y movement automation (PLC based control)
 Kopal fixture

XL - PANEL



1. Program Display
2. Joystick for Y Movement
3. Jog Cycle Start
4. Lamp On/Off
5. Joystick for Z Movement
6. Z-Movement Cycle Start
7. Hydraulic clamp tight/loose
8. Emergency Stop
9. Y Cycle Start

Chennai Metco Pvt Ltd

SP 100A, Ambattur Industrial Estate, Chennai - 600 058, India.
 Ph: 044-2688 0737, 044-2688 0005, Fax: (044) 2688 0850,
 e-mail: mktg@chennaietco.com, web : www.chennaietco.com